

Daylight Savings Time Bulletin

March 5th, 2007

As you may know the dates on which daylight savings time begins and ends will be changing in 2007 as noted below:

- o March 11 (New date for start of daylight savings time)
- o April 1 (Old date for start of DST)
- o October 28 (Old date for end of DST)
- o November 4 (New date for end of DST)

The FDA has posted a notice stating that this change could cause medical devices to behave in an unpredictable way and has asked consumers to request software patches to prevent any problems.

Audioscan equipment is unaffected by the change in daylight savings time.

RM500	Does not contain an internal clock
RM500 SL	Contains an internal clock that does not automatically adjust, the clock can be manually changed by the user.
Verifit	Contains an internal clock that does not automatically adjust, the clock can be manually changed by the user.

If you wish to change the clock on a Verifit or RM500 SL simply select <setup> and pick <date & time> from the menu.

Any questions regarding this issue may be forwarded to:

Dave Cole
Audioscan
41 Byron Avenue
Dorchester ON
N0L 1G0
dave@audioscan.com

audioscan[®]